

Hemoglobin Meter

Mission Hb[®]

*Easy, Fast and Accurate way
for testing of **hemoglobin***

Accurate and Easy to Use

Fast Results (≤ 15 Seconds) of **Hb & HCT**

Capillary & whole blood can be used



IVD

Hb METER

Now Testing of Hemoglobin is easy



- Sample Volume: 10 μ l
- Measuring Time: <15 Seconds
- Battery Life: 1000 Tests
- Calibration: Auto Calibration by Code Chip
- Memory: 1000 Tests Results
- Measuring Range: 5-25 gm/dl

Easy 3 Step Procedure



1. Insert Test Strip



2. Apply Sample



3. Read Result



REMEDY HEALTHCARE PRIVATE LIMITED

Add: S/A-1/20, Jogabai Extn., Near- New Friends Colony, OKHLA, New Delhi-110025. [INDIA].

Ph.: +91-11-65103600, Mb: +91-9911959800/ 9911954800.

Email 1: drcure.poct@gmail.com

Email 2: haemocheck@gmail.com

Web: remedyhealthcare.in

Distributor: