


Gluco Dot[®]

Blood Glucose Monitoring System

No Code System

- High Accuracy
- Easy Operation
- **0.5 μ l** Capillary Blood
- Large LCD
- Biosensor Technology



*5 Seconds
Quick Read*

REMEDY[™]
Good Health For All...

IVD CE 0537

Gluco Dot®

Blood Glucose Monitoring System



Meter Specifications

- **Measuring range:** 20~600 mg/dL
- **Sample Type:** Capillary blood
- **Memory Capacity:** 300 test results with date and time
- **HCT of Tested Blood:** 35%~55%
- **Measurement Unit:** mmol/L or mg/dl
- **Test Time:** 5 Seconds
- **Display:** LCD display
- **Battery Life:** Over 1,000 tests
- **Power supply:** One 3V lithium battery
- **Storage conditions:** Temperature 4°C~30°C
- **Operation Conditions:** Temperature 10°~40°C Humidity Rh10-90%
- **Meter dimensions:** 92x58x19 mm
- **Meter Weight:** 60g

Test Strip Specifications

- **Blood Volume:** 0.5 µL
- **Storage Conditions:** Temperature 4~30°C Humidity Rh≤ 80%

Easy 3 steps
Caring for your health

- Insert test strip to turn on meter



- Apply Blood



- Read Result



5 Seconds
Quick Read

IVD CE 0537

MKT BY:

REMEDY HEALTHCARE PRIVATE LIMITED

Corp. Add: S/A-1/20, Jogabai Extn., Near- New Friends Colony, OKHLA, New Delhi-110025.[INDIA].

Ph.: +91-11-65103600, Mb: +91-9911959800/ 9911954800.

Email 1: remedyhealthcarepvtltd@gmail.com

Web: www.remedyhealthcare.in

REMEDY
Good Health For All...