

REMEDY™

Good Health For All...

(An ISO 9001 : 2015 Certified Company)

HAEMOCHECK™

Haemoglobin Color Scale

“CE”

Registered Product



Certification & Inspection



REMEDY HEALTH CARE PVT. LTD.

(An ISO 9001 : 2015 Certified Company)

Office: S/A-1/20, Jogabai Extn., Near New Friends Colony, Okhla, New Delhi - 110025 (INDIA)

Mob: 8750059800, Help Line No: 9911959800, 9911954800

E-mail: haemochek@gmail.com / remedyhealthcarepvtltd@gmail.com

www.remedyhealthcare.in

WHO REFERENCE

Haemoglobin equivalent g/dl	14	●
	12	●
	10	●
	8	●
	6	●
	4	●

14 g/dl	Healthy
12 g/dl or more	Not Anemic
8-11 g/dl	Mild to Moderate Anemia
6-7 g/dl	Marked Anemia
4-5 g/dl	Severe Anemia
Less than 4 g/dl	Critical

HAEMOCHECK™ : Haemoglobin Color Scale

This Device Can be used in 3 - Simple Steps:



STEP - 1:
Prick the finger with
Sterile blood
Lancet /Needle.



STEP - 2:
Obtain a free flowing
Drop of blood at one end
of Hb-strip and wait for
30 seconds to 1 minute.



STEP - 3:
Match the color of blood
drop from lighter to darker
shade of Color Scale and
note the reading.

- **HAEMOCHECK™** is the simplest device available for estimating the HB level of Anemic Patients in absence of Laboratory under Field Condition in Camps.
- **HAEMOCHECK™** kit contains calibrated, duly validated and “CE” marked Color Shade Card, Specialized HB Test Strips.
- The shades of the Color Scale Card represent the range of HB levels from 4 to 14 g/dl. Just by matching the color of a drop of blood [Fresh Capillary / Venous] on the special test strip provided, you can know the level of Haemoglobin in the blood.
- **HAEMOCHECK™** is best suited for field testing of Haemoglobin of children, pregnant women, Adolescent girl child and people suffering from Malaria or living with other parasitic condition as well as for mass screening in CAMPS.
- Most accurate results than any other method for field testing as per various research & hence Haemoglobin Color Scale Strip based method is recommended by WHO for field testing.
- Haemoglobin Color Scale is being used Internationally as well as by Govt. of India & many state Govt. agencies for Anemia Eradication Programmes viz; National Health Mission, NGO, DHC, PHC, CHC, ANM, ASHA-Workers, Anganwadi-Workers, Medical Practitioners, Clinicians, Gynecologists & Blood Banks etc.
- **HAEMOCHECK™** is Marketed Worldwide by **REMEDY HEALTH CARE PVT. LTD.** An ISO 9001 : 2015 Certified Company.