


Glucose



Specifications		ET-101
		Glucose
Measuring Range		20-600mg/dL (1.1-33.3mmol/L)
Memory Capacity		200 test results
Test Time		10 seconds
Operating Condition	Temperature	57.2-104°F[14-40°C]
	Humidity	≤ 85% humidity
Meter Storage Condition	Temperature	14-140°F[-10-60°C]
	Humidity	≤ 95% humidity
Sample Volume		≥ 4µl
Sample Type		Capillary whole blood
Hematocrit Range		30-55%
Power Supply		1.5V(AAA)x2
Battery Life		More than 1000 times
Display		LCD display (35x45mm)
Dimensions HxWxD(mm)		88 x 64 x 22
Weight		59 grams, without Battery
Technology Used		Amperometric Enzyme Electrode
Time & Date Display		Yes
Glucose Average Result		7 / 14 / 28 Days
Power On & Off		Automatically
Strip Type		Siphonal Strip (Capillary Action)
Warranty		2 Years

Glucose / Hemoglobin



Single Instrument - 2 Parameters

Specifications		ET-222	
		Glucose	Hemoglobin
Measuring Range		20-600mg/dL (1.1-33.3mmol/L)	7-26mg/dL (4.34-16.14mmol/L)
Memory Capacity		200 test results	50 test results
Test Time		10 seconds	10 seconds
Operating Condition	Temperature	57.2-104°F[14-40°C]	
	Humidity	≤ 85% humidity	
Meter Storage Condition	Temperature	14-140°F[-10-60°C]	
	Humidity	≤ 95% humidity	
Sample Volume		≥ 4µl	≥ 8µl
Sample Type		Capillary whole blood	
Hematocrit Range		30-55%	
Power Supply		1.5V(AAA)x2	
Battery Life		More than 1000 times	
Display		LCD display (35x45mm)	
Dimensions HxWxD(mm)		88 x 64 x 22	
Weight		59 grams, without Battery	
Technology Used		Amperometric Enzyme Electrode	
Time & Date Display		Yes	
Glucose Average Result		7 / 14 / 28 Days	
Power On & Off		Automatically	
Strip Type		Siphonal Strip (Capillary Action)	
Warranty		2 Years	


REMEDY™

Healthcare Pvt. Ltd.

Kernel MultiCheck



Meters identify the strips automatically



Glucose Uric Acid
Cholesterol Hemoglobin

Meters and test strips

Mktg. By:

Life's better up here...

REMEDY HEALTHCARE PRIVATE LIMITED

Office : S/A-1/20, Jogabai Extn., Near-New Friends Colony, Okhla, New Delhi-110025. [INDIA].

Ph.No: 011-65103600, Cell: 09911954800.

Email: remedyhealthcarepvtltd@gmail.com

www.remedyhealthcare.in